

# **Non-Targeted Analysis Method Reporting Tool**

## **Instructions for Use**

### **Summary**

The Non-Targeted Analysis (NTA) Method Reporting Tool is a macro-enabled Microsoft Excel workbook that allows for the controlled ontology of method data reporting and the export of the data into a single concise, human-readable file, written in a standard Extensible Markup Language (XML).


### **File Descriptions**

NTA Method Report Tool.xlsm – the base file that can be used to create a new method reporting file.


NTA Method Reporting Tool\_example.xlsm – an example file with completed information.

## Instructions

- 1) Open the NTA Method Report Tool.xlsm, you may have to select “Enable Macros”



- 2) Upon opening, you will be prompted to save the Excel Workbook as a separate file. This file can be saved and edited later prior to exporting the data, so it is recommended to save different files for each individual datafile. Click **Save** to continue.



- 3) The first sheet (titled **Run**) shows the steps for filling out the Method Reporting Tool. This sheet will also tell you if you have completed the minimum amount of information for each section.

|    | A   | B  | C                      | D | E | F | G | H | I | J | K |
|----|---|--|------------------------|---|---|---|---|---|---|---|---|
| 1  | <b>Non-Targeted Analysis Method Reporting Tool</b>  |  |                        |   |   |   |   |   |   |   |   |
| 2  |   |  |                        |   |   |   |   |   |   |   |   |
| 3  | <b>Step 1) Fill out the data in each Input sheet</b>  |  |                        |   |   |   |   |   |   |   |   |
| 4  |   | Sheet  | Complete?              |   |   |   |   |   |   |   |   |
| 5  |   | <a href="#">Sample Input</a>                         | FALSE                  |   |   |   |   |   |   |   |   |
| 6  |   | <a href="#">Chromatography Input</a>                 | FALSE                  |   |   |   |   |   |   |   |   |
| 7  |   | <a href="#">Mass Spectrometry Input</a>              | FALSE                  |   |   |   |   |   |   |   |   |
| 8  |   | <a href="#">QC Method Input</a>                      | FALSE                  |   |   |   |   |   |   |   |   |
| 9  |   | <a href="#">Peaks Input</a>                          | FALSE                  |   |   |   |   |   |   |   |   |
| 10 |   |  |                        |   |   |   |   |   |   |   |   |
| 11 | <b>Step 2) Press "Create Annotation Tables" to annotate the fragments of the peaks [optional]</b> |  |                        |   |   |   |   |   |   |   |   |
| 12 |   | Tables Generated:                                    | FALSE                  |   |   |   |   |   |   |   |   |
| 13 |   |  |                        |   |   |   |   |   |   |   |   |
| 14 | <b>Step 3) Aggregate Data into output file</b>  |  |                        |   |   |   |   |   |   |   |   |
| 15 |   |  | Create Method XML file |   |   |   |   |   |   |   |   |
| 16 |   |  |                        |   |   |   |   |   |   |   |   |
| 17 | <b>Step 4) Give data to NIST</b>  |  |                        |   |   |   |   |   |   |   |   |
| 18 |   | Contact pfas@nist.gov to find out how to share data. |                        |   |   |   |   |   |   |   |   |

- 4) For each section, there is information to fill out under the **Values** heading. For those values that have restricted inputs, there is a drop-down menu.
- a. If there is not an appropriate item in the drop-down menu, contact Ben Place at [benjamin.place@nist.gov](mailto:benjamin.place@nist.gov) with the recommended changes and they will update the form.

| Description    | Value                            |
|----------------|----------------------------------|
| File Name      | PFAC30PAR_PFCA2.mzML             |
| Sample Name    | 3M 3% Lightwater AFFF            |
| Sample Class   | solution reference material      |
| Submitter Name | solution reference material      |
|                | matrix reference material        |
|                | groundwater / surface water      |
|                | drinking water                   |
|                | aqueous film-forming foam (AFFF) |
|                | commercial formulation           |
|                | landfill leachate                |
|                | soil/sediment                    |

- b. When you select a specific value cell, hover text will appear and describe the type of input expected for the value.

| Value  |
|--|
| PFAC30PAR_PFCA2.mzML                           |
| Reference Sta                                  |
| analytical stan                                |
| <a href="mailto:bjp@nist.gov">bjp@nist.gov</a> |

**File Name**  
The full name of the raw data file related to the submitted data. Names should match exactly.

5) The following information is provided for each individual input sheet

**Sheet: Sample\_Input**

| Sample Information |  |  |               |
|--------------------|--|--|---------------|
| User Input         |  |  |               |
| Description        | Value  |  | XML Parameter |
| File Name          | PFAC30PAR_PFCA2.mzML                           |  | filename      |
| Sample Name        | Reference Standard for PFAS                    |  | description   |
| Sample Class       | analytical standard                            |  | class         |
| Submitted by       | <a href="mailto:bjp@nist.gov">bjp@nist.gov</a> |  | submitter     |

*File Name:* The name of the raw data file to be submitted, can be an mzML or proprietary data format. Names should match exactly.

*Sample Name:* name of the sample as described by the user, uncontrolled.

*Sample Class:* the category or class to which the sample belongs.

Examples: analytical standard, landfill leachate

*Submitter Name:* Unique identifier for submitting user or laboratory. For individual users, use email address.

Sheet: Chromatography\_Input

|    | A  | B                       | C | D                    |
|----|--|-------------------------|---|----------------------|
| 1  | <b>Chromatography</b>                              |                         |   |                      |
| 2  | <b>User Input</b>                                  |                         |   |                      |
| 3  | <b>Description</b>                                 | <b>Value</b>            |   | <b>XML Parameter</b> |
| 4  | Type of Chromatography                             | Liquid Chromatography   |   | ctype                |
| 5  | Chromatograph Vendor                               | ThermoFisher Scientific |   | cvendor              |
| 6  | Chromatograph Model                                | UltiMate 3000           |   | cmodel               |
| 7  | Sample Solvent                                     | water                   |   | ssolvent             |
| 8  | Mobile Phase 1 Solvent                             | water                   |   | mp1solvent           |
| 9  | Mobile Phase 1 Additive                            | ammonium acetate        |   | mp1add               |
| 10 | Mobile Phase 2 Solvent                             | water                   |   | m2solvent            |
| 11 | Mobile Phase 2 Additive                            | ammonium acetate        |   | mp2add               |
| 12 | Mobile Phase 3 Solvent                             | none                    |   | mp3solvent           |
| 13 | Mobile Phase 3 Additive                            | none                    |   | mp3add               |
| 14 | Mobile Phase 4 Solvent                             | none                    |   | mp4solvent           |
| 15 | Mobile Phase 4 Additive                            | none                    |   | mp4add               |
| 16 | Analytical Column Vendor                           | Agilent Technologies    |   | colvendor            |
| 17 | Analytical Column Brand Name                       | Poroshell C18           |   | colname              |
| 18 | Analytical Column Phase                            | C18                     |   | colchemistry         |
| 19 | Analytical Column Inner Diameter                   | 2.1                     |   | colid                |
| 20 | Analytical Column Length                           | 150                     |   | collen               |
| 21 | Analytical Column Particle Diameter/Film Thickness | 1.8                     |   | coldp                |
| 22 | Guard Column Vendor                                | none                    |   | gcolvendor           |
| 23 | Guard Column Brand Name                            |                         |   | gcolname             |
| 24 | Guard Column Phase                                 | none                    |   | gcolchemistry        |
| 25 | Guard Column Inner Diameter                        |                         |   | gcolid               |
| 26 | Guard Column Length                                |                         |   | gcollen              |
| 27 | Guard Column Particle Diameter/Film Thickness      |                         |   | gcoldp               |
| 28 |  |                         |   |                      |

*Type of Chromatography:* the type of chromatography used for the generation of the data.

Options: liquid chromatography (LC), gas chromatography (GC), capillary electrophoresis (CE), or none

*Chromatograph Vendor:* the manufacturer or vendor of the chromatograph

Options: most common chromatograph vendors

*Chromatograph Model:* the model name of the chromatograph

Example: UltiMate 3000, Agilent 1260, Waters 2695

*Sample Solvent:* the primary solvent of the sample, as injected

Options: most common solvents

Example: if sample solvent is 30 % water, 69 % methanol, and 1 % formic acid, the primary solvent is methanol

*Mobile Phase Solvent 1/2/3/4:* the primary solvent used for the mobile phase program, allowing for up to a quaternary pump. Select 'none' if not used in program.

Options: most common solvents.

Example: if Mobile Phase 1 (A) is 95 % water, 4.9 % acetonitrile, and 0.1 % formic acid, the primary solvent is water

*Mobile Phase Additive 1/2/3/4:* the primary additive used in the respective solvent for the mobile phase program, allowing for up to a quaternary pump. Select 'none' if not used in program.

Options: most common additives

Example: if Mobile Phase 2 (B) is 95 % acetonitrile with water and 10 mM ammonium acetate, adjusted to pH 6 with ammonium hydroxide, the primary additive is ammonium acetate.

*Analytical/Guard Column Vendor:* the manufacturer of the analytical or guard chromatography column used for the experiment.

Options: most common column manufacturers

*Analytical/Guard Column Brand Name:* the commercial brand/product name of the analytical or guard chromatography column.

Examples: Halo C18, Eclipse Plus C18

*Analytical/Guard Column Phase:* the ligand chemistry of the stationary phase in the analytical or guard column.

Options: C18, C8, pentafluorophenyl (PFP), biphenyl, silica, mixed-phase, diol, pentadiol

Note: if the phase has more than one chemistry or multiple columns in sequence are used, select 'mixed-phase'

*Analytical/Guard Column Inner Diameter/Length:* the dimensions of the analytical or guard column.

Note: for LC columns, report values in millimeters. For GC or CE columns, report values in micrometers.

*Analytical/Guard Column Particle Diameter/Film Thickness:* the particle diameter (LC) or the film thickness (GC, CE) of the analytical or guard column.

Note: for LC columns, report particle diameter in micrometers. For GC or CE columns, report film thickness in micrometers.

Sheet: Mass Spectrometry\_Input

|    | A                                | B  | C | D                    |
|----|----------------------------------|--|---|----------------------|
| 1  | <b>Mass Spectrometry</b>         |  |   |                      |
| 2  | <b>User Input</b>                |  |   |                      |
| 3  | <b>Description</b>               | <b>Value</b>                             |   | <b>XML Parameter</b> |
| 4  | Mass Spectrometer Vendor         | ThermoFisher Scientific                  |   | msvendor             |
| 5  | Mass Spectrometer Model          | Q-Exactive                               |   | msmodel              |
| 6  | Ionization Mode                  | atmospheric pressure chemical ionization |   | imode                |
| 7  | Polarity                         | negative                                 |   | polarity             |
| 8  | Ionization Voltage/Current       | 4500                                     |   | vvalue               |
| 9  | Ionization Voltage/Current Units | V  |   | vunits               |
| 10 | First Mass Analyzer              | quadrupole                               |   | massanalyzer1        |
| 11 | Second Mass Analyzer             | quadrupole                               |   | massanalyzer2        |
| 12 | Fragmentation Type               | HCD                                      |   | fragmode             |
| 13 | Fragmentation Energy             | 15;30;45                                 |   | cevalue              |
| 14 | Fragmentation Energy Type        | stepped                                  |   | cetype               |
| 15 | Fragmentation Energy Units       | normalized                               |   | ceunits              |
| 16 | MS2 Experiment                   | DDA                                      |   | ms2exp               |
| 17 | Isolation Width/Window (Da)      | 0.5                                      |   | isowidth             |
| 18 | Instrument Mass Accuracy (ppm)   | 5  |   | msaccuracy           |
| 19 | MS1 Resolution                   | 30000                                    |   | ms1resolution        |
| 20 | MS2 Resolution                   | 15000                                    |   | ms2resolution        |

*Mass Spectrometer Vendor:* the manufacturer or vendor of the mass spectrometer

Options: most common mass spectrometer vendors

*Chromatograph Model:* the model name of the chromatograph

Example: Q-Exactive HF, Agilent 6550, Sciex X500

*Ionization Mode:* the type of ionization mechanism used with the mass spectrometer.

Options: atmospheric pressure chemical ionization, atmospheric pressure photoionization, electrospray ionization, electron ionization, chemical ionization

*Polarity:* the ionization polarity used for the compounds reported.

Options: positive, negative

Note: Only one polarity can be selected at a time. If polarity switched was used, the worksheet must be submitted twice (once for positive, once for negative).

*Ionization Voltage/Current:* the voltage or current applied for the ionization mode.

Note: Report number only, units will be selected in the next section.

*Ionization Voltage/Current Units:* the units for the reported value of ionization voltage/current.

Example: V,  $\mu$ A



*First Mass Analyzer:* the first mass analyzer of the mass spectrometer. If only one mass analyzer is used, then this is the only mass analyzer used.

Options: quadrupole, linear ion-trap, magnetic sector, time-of-flight, orbitrap, fourier transform ion cyclotron resonance.

Example: for a QTOF instrument, the first mass analyzer is a quadrupole.

*Second Mass Analyzer:* the second mass analyzer of the mass spectrometer. If only one mass analyzer is used, then select 'none'.

Options: quadrupole, linear ion-trap, magnetic sector, time-of-flight, orbitrap, fourier transform ion cyclotron resonance.

Example: for a LIT-Orbitrap, the second mass analyzer is an orbitrap.

*Fragmentation Type:* for the generation of fragment ions, report the type of fragmentation or collision cell used.

Options: high-energy collisional dissociation (HCD), collisionally-induced dissociation (CID), electron capture dissociation (ECD), in-source fragmentation

Note: For data using GC-EI-MS, select 'in-source fragmentation'.

*Fragmentation Energy:* the energy value used for fragmentation. For stepped/ramped collision energy, separate levels by a semicolon (;). The type of fragmentation energy and units will be selected below.

*Fragmentation Energy Type:* the type of fragmentation energy used. For stepped/ramped collision energy, selected 'stepped'. For all other types of fragmentation energy, select 'fixed'.

*Fragmentation Energy Units:* the units of the reported energy value used for fragmentation.

Options: normalized, volts (V)

*MS2 Experiment:* the type of fragmentation (MS2) experiment used for generating fragment ions.

Options:

DDA: For data-dependent analysis, where a narrow window for a specific precursor ion is selected, including TopN and Information Dependent Acquisition (IDA).

SWATH: For sequential windowed-based precursor selection, including the Sequential Windows of All Theoretical Mass Spectra (SWATH) technique

DIA: For data-independent analysis, where no precursor ion is selected, including All Ion Fragmentation or MS<sup>E</sup> analysis.

none: if using in-source fragmentation

*Isolation Width/Window:* The isolation width or window for the precursor ion selection.

Note: for DDA, report the isolation width in Da. For SWATH, report the window width in Da. For DIA or none, leave blank.

*Instrument Mass Accuracy:* the mass accuracy of the instrument, reported in ppm.

Note: if you are not reporting data using high resolution mass spectrometry, enter 0.

*MS1 Resolution:* the resolution of the MS1 data.

*MS2 Resolution:* the resolution of the MS2 data. If no MS2 data, leave empty.

Sheet: QC Method\_Input

|   | A                              | B            | C | D                    |
|---|--------------------------------|--------------|---|----------------------|
| 1 | <b>QC Method</b>               |              |   |                      |
| 2 | <b>User Input</b>              |              |   |                      |
| 3 | <b>Description</b>             | <b>Value</b> |   | <b>XML Parameter</b> |
| 4 | Was a QA/QC Method Used        | TRUE         |   | qcused               |
| 5 | Mass Analyzer Calibration      | TRUE         |   | qctype               |
| 6 | External Standard Verification | TRUE         |   | qctype               |
| 7 | Internal Standard Verification | FALSE        |   | qctype               |
| 8 | Matrix Standard Verification   | FALSE        |   | qctype               |

*Was a QA/QC Method Used?* If the user has established (internally or published) protocols for the quality assurance and/or quality control for the measurement of the submitted data, select TRUE.

Note: if your QA/QC protocol is not described in the below techniques, but you still used a QA/QC protocol, select TRUE.

*Mass Analyzer Calibration:* Was the mass analyzer used calibrated within the period specified by the manufacturer? If yes, select TRUE.

*External Standard Verification:* Was the quality of the instrumental method verified by the analysis of a solvent spiked with known chemical standards during the analytical sequence that produced the submitted data? If yes, select TRUE.

Example: an aliquot of methanol was spiked with perfluoroalkyl carboxylic acids and analyzed during the same sequence as the samples.

*Internal Standard Verification:* Was the quality of the instrumental method verified by the analysis of samples that have been enriched/spiked with known chemical standards? If yes, select TRUE.

Example: the sample was spiked with <sup>13</sup>C-labeled perfluoroalkyl carboxylic acid standards prior to being analyzed.

*Matrix Standard Verification:* Was the quality of the instrumental method verified by the analysis of a control sample, which is a material consisting of a matrix similar to the unknown samples that contains known chemical compounds that are endogenous or enriched/spiked? If yes, select TRUE.

Example: a certified reference material of fish tissue extracted into methanol and analyzed during the same sequence as the samples.

## Sheet: Peaks\_Input

|    | A                         | B                                      | C                | D         | E         | F                     | G                   | H                |
|----|---------------------------|--|------------------|-----------|-----------|-----------------------|---------------------|------------------|
| 1  | <b>Mass Spectrometry</b>  |  |                  |           |           |                       |                     |                  |
| 2  | <b>User Input</b>         |  |                  |           |           |                       |                     |                  |
| 3  | <b>name</b>               | <b><a href="#">Suspect List ID</a></b> | <b>Ion State</b> | <b>mz</b> | <b>rt</b> | <b>peak_starttime</b> | <b>peak_endtime</b> | <b>Verified?</b> |
| 4  | Perfluorobutanoic acid    | 2649                                   | [M-H]-           | 212.9792  | 4.5       | 4.1                   | 4.8                 | TRUE             |
| 5  | Perfluoropentanoic acid   | 2646                                   | [M-H]-           | 262.976   | 8.6       | 8.4                   | 8.9                 | FALSE            |
| 6  | Perfluorohexanoic acid    | 2643                                   | [M-H]-           | 312.9728  | 10.8      | 10.5                  | 11.1                | FALSE            |
| 7  | Perfluoroheptanoic acid   | 2640                                   | [M-H]-           | 362.9696  | 12.09     | 11.9                  | 12.35               | FALSE            |
| 8  | Perfluorooctanoic acid    | 2637                                   | [M-H]-           | 412.9664  | 13.05     | 12.8                  | 13.3                | TRUE             |
| 9  | Perfluorononanoic acid    | 2635                                   | [M-H]-           | 462.9632  | 13.8      | 13.6                  | 14.1                | FALSE            |
| 10 | Perfluorodecanoic acid    | 2632                                   | [M-H]-           | 512.96    | 14.5      | 14.3                  | 14.8                | FALSE            |
| 11 | Perfluorododecanoic acid  | 2629                                   | [M-H]-           | 612.9537  | 15.6      | 15.4                  | 15.8                | TRUE             |
| 12 | Perfluoroundecanoic acid  | 2630                                   | [M-H]-           | 562.9568  | 15.1      | 14.9                  | 15.4                | TRUE             |
| 13 | Perfluorotridecanoic acid | 2628                                   | [M-H]-           | 662.9505  | 16        | 15.8                  | 16.3                | TRUE             |

To submit mass spectra for identified compounds within a single raw data file, you must provide:

*Name:* name of the compound. This is an uncontrolled value and does not need to match the compound name in the NIST Suspect List.

*Suspect List ID:* the ID value for the compound in the NIST Suspect List of Possible Per- and Polyfluoroalkyl Substances (<https://data.nist.gov/od/id/mds2-2387>)


*Ion State:* the ion/adduct state of the precursor ion as it related to the measured  $m/z$  value for the specific compound.

*mz:* the measured  $m/z$  for the precursor ion used for the fragmentation of the compound, not the exact mass of the compound. For in-source fragmentation, this is  $m/z$  corresponding the molecular or pseudomolecular ion of the compound.

*rt:* the retention time of the chromatographic peak apex for the identified compound. This can be reported from most software, or can be approximated by the user.

*peak\_starttime:* the retention time of the start of the chromatographic peak for the identified compound. This can be reported from most software, or can be approximated by the user.

*peak\_endtime:* the retention time of the end of the chromatographic peak for the identified compound. This can be reported from most software, or can be approximated by the user.




*Verified?* If the compound identity was verified using a chemical standard, by matching the retention time, precursor  $m/z$  and/or fragmentation mass spectrum, select TRUE.

Note: if the sample is a chemical standard for the identified compound, select TRUE.

- 6) When the information is complete, return to the **Run** sheet.
  - a. If you have peaks for which you want to annotate the fragments, go to step 7.
  - b. If you have no peaks to annotate, go to step 9.

- 7) To annotate fragments, press the **Create Annotation Tables** button.



For each compound identified in the **Peaks\_Input** sheet, a new sheet will be generated titled “Ann\_[compound name]”.

- 8) For each compound, you will be able to submit annotated fragments for the fragmentation mass spectrum. Annotation is the attribution of the elemental formula (and structure, if possible) to a specific measured  $m/z$  of a fragment. All fragments within a mass spectrum do not need to be annotated. All compounds submitted with the dataset do not need to have annotated fragments.

|   | A   | B                          | C                              | D        | E                           | F |
|---|---|----------------------------|--------------------------------|----------|-----------------------------|---|
| 1 | <b>Peak Annotation for Perfluorohexanoic acid</b> |                            |                                |          |                             |   |
| 2 | User Input  |                            |                                |          |                             |   |
| 3 | Fragment $m/z$                                    | Fragment Elemental Formula | Fragment SMILES                | Radical? | Fragment Citation           |   |
| 4 | 118.9926 C2F5                                     |                            | F[C-](F)C(F)(F)F               | FALSE    | DOI:10.1002/rcm.3274        |   |
| 5 | 268.983 C5F11                                     |                            | FC(F)(C(F)(F)C(F)(F)C(F)(F)[C- | FALSE    | DOI:10.1002/rcm.3274        |   |
| 6 | 59.0139 C2H3O2                                    |                            | [O-]\\C=C\\O                   | FALSE    | DOI:10.1021/acs.est.6b05843 |   |
| 7 |   |                            |                                |          |                             |   |
| 8 |   |                            |                                |          |                             |   |
| 9 |   |                            |                                |          |                             |   |

For each fragment for each compound, you can submit the following:

*Fragment  $m/z$*  – the measured  $m/z$  value for the fragment to-be-annotated (required)

*Fragment Elemental Formula* – the elemental formula attributed to the measured fragment  $m/z$  (required)

*Fragment SMILES* – the chemical structure (in SMILES notation) attributed to the measured fragment  $m/z$  (optional)

*Radical?* If the annotated fragment contains a radical electron, enter TRUE. If not, enter FALSE. If the existence of a radical electron cannot be determined, enter UNKNOWN.

*Fragment Citation* – the documented evidence of the elemental formula, structure, and/or radical for the measured fragment  $m/z$ .

Options:

DOI – if the evidence is listed in a published manuscript, enter the DOI with the format:  
DOI:number

Website: if the evidence is listed on a website, enter the URL with the format:  
<https://www.website.com/12345>

User interpretation: if the user interpreted the elemental formula, structure, or radical without supporting evidence, enter USER.

Note: If software was used to annotate the fragments, enter the software website or DOI as the citation. For example: <https://cfmid.wishartlab.com/>

- 9) When completed, check to make sure all parts of Step 1 are completed by going to the **Run** sheet.

|    |  |   |                  |  |
|----|--|---|------------------|--|
| 1  | <b>Method Reporting Tool for Non-Targeted Analyses</b> |   |                  |  |
| 2  |  |   |                  |  |
| 3  | <b>Step 1) Fill out the data in each Input sheet</b>   |   |                  |  |
| 4  |  | <b>Sheet</b>                            | <b>Complete?</b> |  |
| 5  |  | <a href="#">Sample Input</a>            | TRUE             |  |
| 6  |  | <a href="#">Chromatography Input</a>    | TRUE             |  |
| 7  |  | <a href="#">Mass Spectrometry Input</a> | TRUE             |  |
| 8  |  | <a href="#">QC Method Input</a>         | TRUE             |  |
| 9  |  | <a href="#">Peaks Input</a>             | TRUE             |  |
| 10 |  |   |                  |  |

If all sheets have TRUE under complete, you can press **Create Method XML File**

|  |                               |
|--|-------------------------------|
| <b>Step 3) Aggregate Data into output file</b> | <b>Create Method XML file</b> |
|  |                               |

If successful, a pop-up message should notify you of the name of the new XML file, which will be located in the same folder as the NTA Method Reporting Tool.xlsm file.

